

## 5L Oxygen Concentrator

Model: HB512

### Main Features



Timing accumulates function



Durable High-Effective molecular sieve



4 alarm systems:  
-low oxygen concentration alarm  
-power off alarm  
-pressure circulation fault alarm  
-compressor fault alarm



Continuous operation for 24 hour/day.

### Suitable range:

Clinic, homecare, nursing house, geriatric nursing etc

### For the crowd:

Geriatrics, pregnant women, high mental labor, students, white-collar, and peoples suffer chest distress, lapse of memory, absent-minded, indigestion, bad sleeping etc.

### Specifications

Model No	HB512
Oxygen flow	1—5 LPM
Oxygen purity	93% ± 3%
Low oxygen alarm	No
Optional	Nebulizer, 2-way splitter
Power consumption	<320 W
Output pressure	0.05–0.07Mpa
Operation noise(dB)	< 45dB ( A )
Net weight	35 lbs (15.9 kg)
Operating temp	5°C —40°C
Dimensions (L x W x H)	43 x 32 x 62 cm
Equip. Classification	Class II Type B
Power Supply	AC 220V ± 22V/50Hz or 110V ± 10V /60Hz



### Brief Introduction:

They require very few maintenance among their lifetime. All inner components are modularization designed, very easy to be maintained by general persons. The work performance are very steady.